3D Guidance drive BAY™

Class 1, Type B Applied Part



Track Objects with Magnetic DC Technology

- ➤ Fast, dynamic tracking as many as 420 updates per second.
- Miniaturized passive sensors outputs immune to "power-line" noise sources.
- Compact packaging electronics conveniently fit into drive bay of your PC.
- ➤ All attitude tracking no inertial drift or optical interference.
- High metal immunity no distortion from non-magnetic metals.







Magnetic field transmitter options for mid and shortrange tracking

3D Guidance drive BAY™

Technical

Sensor Configurations

Degrees of Freedom **Update Rate**

Translation Range

Angular Range

Static Accuracy*

Static Resolution

Model 55 (0.56 mm), Model 90 (0.9 mm), Model 130 (1.5 mm), Model 180 (2.0 mm), Model 800 (8.0 mm)

6 (Position and Orientation)

Up to 420 updates/second for each sensor (Default: 240 updates/second)

MODEL 55 SENSOR

- Mid-Range Transmitter: 25 cm (10.0 inches)
- Short-Range Transmitter: Contact Ascension MODEL 90 SENSOR
- Mid-Range Transmitter: 36 cm (14.0 inches)
- Short-Range Transmitter: Contact Ascension **MODEL 130 SENSOR**
- Mid-Range Transmitter: 46 cm (18.0 inches)
- Short-Range Transmitter: Contact Ascension **MODEL 180 SENSOR**
- Mid-Range Transmitter: 58 cm (23.0 inches)
- Short-Range Transmitter: Contact Ascension
- **MODEL 800 SENSOR** • Mid-Range Transmitter: 78 cm (31.0 inches)
- Short-Range Transmitter: 46 cm (18.0 inches)

All Attitude: ± 180° Azimuth & Roll, ± 90° Elevation Position: 1.4 mm (0.055 inch) RMS Orientation: 0.5° RMS

*Higher accuracies achievable in smaller tracking volumes

Position: 0.5 mm (0.02 inch) at 30.5 cm (12.0 inches) Orientation: 0.1° at 30.5 cm (12.0 inches) *Resolution measured for tracker with mid-range transmitter and 8 mm sensor.

X, Y, Z positional coordinates, orientation angles, orientation matrix or quaternions

USB 2.0

Binary data records Windows API and Drivers

Physical

Data Format

Communication

Outputs

Interface

Electronics Unit

Transmitters

Passive Sensors

Power

Environment

Operating Temperature

20.0 cm (7.8 inches) x 14.7 cm (5.8 inches) x 4.1 cm (1.6 inches) metal box fitting into a 13.3 cm (5.25 inches) PC drive bay

- Mid-Range: 9.6 cm (3.8 inches) cube with 3.3 m (10.9 ft.) cable
- Short-Range: 5.3 cm (2.9 inches) x 5.3 cm (2.9 inches) x 6.9 cm (2.7 inches) with 3.3 m (9.8 ft.) cable

MODEL 55: 0.56 mm (0.02 inch) x 80 mm (3.2 inches) to 210 mm (8.27 inches) for biopsy needle configurations only, with 2.3 m (7.5 ft.) cable

MODEL 90: 0.9 mm (0.04 inch) x 7.25 mm (0.29 inch) with 3.3 m (10.9 ft.) cable

MODEL 130: 1.5 mm (0.05 inch) x 7.7 mm (0.30 inch) with 3.3 m (10.9 ft.) cable

MODEL 180: 2.0 mm (0.07 inch) x 9.9 mm (0.38 inch) with 3.3 m (10.9 ft.) cable

MODEL 800: 8.0 mm (.31 inch) x 20.0 mm (0.78 inch) with 3.3 m (10.9 ft.) cable

- Model 55, 90, 130 & 180 only:
 Ascension Medi-Mag Cable, USP class 6 jacket material.
- USP class 6 sensor housing.

Assembly and cable materials are EtO and cold sterilant tolerant. Warning: Semiconductor devices in sensor connector are not gamma shielded and may be damaged or erased if exposed to gamma radiation and/or autoclaving.

Sensors and cable assemblies are fragile components and must be sheathed, isolated and safeguarded prior to use in patients.

Drawn from PC's power supply through rear panel connector.

- + 12V: 1.6A nominal; 2.9A maximum
- + 5V: 600mA nominal

5°C to 40°C; 90% non-condensing humidity

Ferromagnetic objects and stray magnetic fields in the operation volume may degrade performance. Contact us for assistance in minimizing metallic distortion and noise interference.

FEATURE BENEFITS

Metal tolerant

Outputs unaffected by composite materials. Capable of driving errors induced by highly conductive metals (such as aluminum) to zero by adjusting measurement rate.

Advanced new magnetic technology and signal processing • Improved dynamic performance over longer ranges.

• "Power-line" noise filtered out.

Occlusion and drift free Clear line-of-sight between transmitter and sensor(s) is not required.

Body mountable transmitter New lightweight coil set can be mounted on head or body.

Onboard diagnostics

Self-diagnostics and run-time monitoring for improved tracker reliability and safety.

Software support

USB tracker control API for XP/Pro, XP, Vista, Window 7, 32 & 64 bit, Sample

programs.



Tracker electronics fit into the drive bay of a PC chassis and uses its power. No additional power supply required.



Regulatory Certifications

- Class I Medical Device with Type B Applied Part (Sensors).
- RoHS and WEEE compliant.
- Medical users must comply with all pertinent FDA/CE/IRB certifications prior to using this device in human patients.

Note on Accuracy

Accuracy is defined as the root mean square (RMS) deviation of a true measurement of the magnetic center of a single sensor with respect to the magnetic center of a single transmitter measured over the specified translation range. Accuracy varies from one location to another over this range and will be degraded if there are interfering electromagnetic noise sources or metal in the operating environment, which have not been identified and minimized



© 2012 Ascension Technology Corporation / 3D Guidance driveBAY is an Ascension Trademark 11/12